
Instruction for use



Ambu® Man Advanced

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
www.ambu.com

492 2844 40 - 2016/01 - V01 Printed in China

English Directions for use	日本語 使用法
1	121
Česky Navod k použití	Nederlands Gebruiksaanwijzing
13	133
Dansk Brugsanvisning.....	Norsk Brukerveiledning
25	145
Deutsch Bedienungsanleitung	Polski Instrukcja obsługi
37	157
Ελληνικά δηγίες Χρήσεως	Portugues Manual de instrucoes
49	169
Espanol Manual de instrucciones	Русский Directions for use
61	181
Suomi Kayttoohje.....	Svenska Instruktionshandbok.....
73	195
Francais Mode d’emploi	Turkce Kullanım talimatları
85	207
Magyar Hasznalati utmutato.....	中文 使用指南
97	219
Italiano Manuale d’uso	
109	

Safety and Regulatory Notices

FCC Statement

This device complies with part 15 of the FCC Rules, Operation is subject to the following two conditions:

- 1) This device may not cause harmful interference and,
- 2) This device must accept any interference received, including interference that may cause undesired operations

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications to the product are not allowed.

Contains FCC ID: **PD98260NG**

Content	Page
1. Introduction	3
2. Restrictions and Cautions	3
3. Specifications	4
4. Functions	5
5. Preparation for training	7
6. Usage of the manikin	8
7. Cleaning	10
8. Quick connection guide	11

Ambu® is a registered trademark of Ambu A/S, Denmark.
Ambu is certified according to ISO 9001 and ISO 13485.

This product complies with the essential requirements of Directive 1999/5/EC of the European Parliament and of the Council of 9 March 1999 on radio equipment and telecommunications terminal equipment and the mutual recognition of their conformity.

1. Introduction

The Ambu® Man Advanced is an advanced instruction and training manikin for the simulation of realistic conditions during cardiopulmonary resuscitation by means of securing the airway giving chest compression and performing defibrillation. Moreover the AmbuMan Advanced will cover all aspect of advanced life-support training.

The system is designed for practising the implementation and co-ordination of vital stages of the resuscitation process using a advanced training manikin that can simulate the most important ECG signals for advanced life-support training.

The manikin has the size of a normal adult and provides an exceptionally lifelike representation of the human anatomy, particularly of those features important to training in modern resuscitation techniques.

The built-in instrumentation shows the effectiveness of the resuscitation, including results achieved for insufflation volume and depth of external chest compression, and indicates any stomach inflation and incorrect hand positioning as well.

2. Restrictions and Cautions

Wideband Data Transmission Systems

2400,0 – 2483,5 MHz

The product can be used in EU member states and EFTA countries in respect to the following restrictions.

Country	Restriction
France	Not allowed for outdoor use.
Italy	Not allowed for outdoor use.
Luxembourg	Implemented general authorisation is required for network and service supply.
Norway	Not allowed to be used within a radius of 20 km from the centre of Ny-Ålesund.
Russian Federation	Not allowed to be used before national approval based on the national standard system (GOST) and conformity certificate.

Caution

Only use the manikin in dry surroundings. Do not expose the manikin directly to any kind of liquid.

Remove the battery pack if the manikin is not used for a longer period of time.

If you detect any smoke or smell from the manikin turn off the manikin immediately and stop using the manikin.

If the manikin has been stored in the cold, let the manikin warm to room temperature before use and assure that no condensation has occurred as this could harm the electronics.

Do not interchange the defibrillation electrodes and ECG electrodes on the chest of the manikin.

When using a defibrillator, it is important to follow the manufacturer's instructions for use, particularly the safety rules.

- Before starting the actual training, the instructor should inform the trainees about the proper use of the defibrillator, alerting them to the hazards involved in disregarding the safety rules.
- The manikin is sensitive to moisture. Make sure it does not get wet, as this might cause a short-circuit.

3. Specifications

3.1. Weight:

Full body with carrying cases: approx. 13 kg

Full body in portable suitcase: approx. 26 kg

3.2. Dimension:

Full body ready for training: approx. 170 cm

Full body in portable suitcase (L x W x H): approx. 105 x 62 x 25 cm

3.3. Part/Material

Basic unit	Polyethylene
Instrument part	ABS plastic
Skull	PVC, hard
Chest skin	PVC, soft
Face piece	PVC, soft
Training suit	50% cotton and 50% polyester

3.4. Battery

The AmbuMan Wireless can be equipped with an battery pack.

3.5. Battery lifetime

The battery lifetime using WLAN is approximately 10 hours using fully charged battery pack.

3.6. Adapter (optional)

AC/DC Adapter

Input: 100-240 V AC / 47-63 Hz / 700 mA)

Output: 12 V DC / 2,0 A

3.7. Operating distance

The wireless connection has a range of approximately 50 meters outside having no obstacles in between.

3.8. Usage and storage:

Storage temperature (non-condensing): -18 °C (-0,4 °F) to 40 °C (105 °F)

Usage temperature (non-condensing): -5 °C (23 °F) to 40 °C (105 °F)

Humidity: 5% to 95%

4. Functions

4.1. Defibrillation and ECG ①

The manikin is equipped with electronics that can simulate an ECG-Signal on the manikin. The control of the signal is done by the virtual ECG-Box inside the Ambu CPR-Software or by using the external Ambu ECG-Box.

4.1.1. Defibrillation electrodes

The manikin is equipped with two defibrillation electrodes in order defibrillate the manikin. The defibrillation electrodes (1.1) are placed on the chest of the manikin on the right chest next to the upper sternum and on left chest around the 5th to 6th rip. (1.2)

If semi-automatic defibrillators or AED's with adhesive chest electrodes are used, special adapters for the specific defibrillator (simulating the adhesive chest electrodes) should be fitted to the chest defibrillation electrodes on the manikin making it possible to connect the cable from the AED's.

4.1.2. ECG electordes

4 ECG Electrodes (1.3) are placed on the chest of the manikin making it possible to monitor the ECG on a defibrillator or a separate scope.

4.2. Monitoring instrument ②

To pull out the monitoring instrument depress the catch (2.1) and let the instrument slide out. The training can also be performed with the monitoring instrument inside the manikin.

Readings on the monitoring instrument can be viewed from two sides. By activating the cover plate on the side facing the trainee, instrument readings can be concealed from trainees undergoing a CPR test (2.2), while allowing the training instructor to monitor the effectiveness of the CPR on the opposite side.

The monitoring instrument is divided into two halves. The left-hand side registers ventilation, including inspiration volume (2.3), and provides indication of stomach inflation (2.4). The right-hand side of the instrument registers chest compression, including depth of compression in millimetres (2.5), and provides indication of incorrect hand positioning (2.6).

The gauges for insufflated volume and depth of compression are designed to display green when the correct volume and chest compression have been reached, and red when the measurement recorded is outside the correct range.

If the display is green and red the measurement is at the borderline of the correct range. The instrument registers inflation in the stomach (2.4) and incorrect hand positioning (2.6) by changing from black to red colour display.

The values on the scale plates on the monitoring instrument are in accordance with the actual guidelines for resuscitation. It is possible to order instrument plates with other recommendations or guidelines.

4.2.1. LAN connector ②



To connect the manikin to a computer the LAN connector (2.7) can be used, placed on backside of the manikin.

4.2.2. Power button ②

To turn on the manikin for the use of the software application press the power button (2.8) on top of the monitoring instrument. The green LED (2.9) is flashing while the system is booting and is steady activated once the system is completely powered up.

To turn off the manikin the power button needs to be pressed at least 3 seconds. If the button is pressed for 10 seconds and more the system is forced to shut down.

An overview of the different LED status can be seen in the following table:

LED status	LED	Description	Flashing pattern
Green steady	●	System running	
Green slow flashing	○	System booting	
Off	○	System off	

4.2.3. Reset to default settings

To reset the manikin into the default settings press reset button (2.10) using a pin.

5. Preparation for training

The manikin is supplied in a portable suitcase. Open the suitcase and take out the manikins and legs.

5.1. Mounting the legs on the torso ③

Pull down the trousers a little around the hips. Place the legs in a position so that the 2 dowels on the hip part can slide into the corresponding dents at the bottom of the manikin's body, see 3.1.

Press the two Velcro straps firmly to the body, see 3.2.

To remove legs pull off the Velcro straps and the legs will be disengaged.

5.2. Monitoring Instrument

Activate the instrument by depressing the catch. In case the instrument is not activated (pulled out), no damage is done to the manikin nor to the instrument during chest compression.

5.3. Adjusting chest stiffness ④

The stiffness of the chest can be adjusted as required by loosening the thumb screw on the back of the manikin: For reduced stiffness set to 'LOW' position, for increased stiffness set to 'HIGH' position.

The values shown, approximately 6 N/mm (0,6 kg/mm) and 11 N/mm (1,1 kg/mm), indicate the force, which must be applied to compress the chest by 1 mm.

Example: Compressing the chest by 40 mm at the 'LOW' setting, a force of approximately 240 Newton (24 kg) must be applied.

The normal setting is 'MEDIUM' corresponding to approximately 8,5 N/mm (0,85 kg/mm).

5.4. Battery placement ⑤

The battery compartment is placed on the backside of the manikin. Depress the catch to open (5.1) and remove the battery pack.

NOTE: Take care that the battery pack is placed correctly in the battery compartment.

5.5. Power supply ⑤

The manikin can be powered using a universal external AC/DC adapter (output 12 V / 2,0 A) with a straight connection plug (5,5x2,1x12 mm). The connector can be connected to the socket place inside the battery compartment, see 5.2.

6. Usage of the manikin

6.1. Ventilation ⑥

The manikin is not supplied with a hygienic system, as normally a resuscitator is used. Mouth-to-mouth or mouth-to-nose ventilation must therefore not be performed.

Ventilation can be performed by using a resuscitator with mask (6.1), ventilation through an endotracheal tube (6.2) or other modern airway management devices.

Correct ventilation makes the manikin's chest raise and fall. When doing bag/mask ventilation there is a possibility for stomach inflation, which causes the manikin's stomach to rise. This happens if the ventilation is too fast, the volume is too high or if the head is not tilted correctly. As in real-life situations ventilating through the intubation tube eliminates the risk of stomach inflation if the tube is correctly placed.

6.2. Intubation

The manikin can be intubated through the mouth (oral intubation). Intubation through the nose (nasal intubation) is not possible.

The intubation can be performed with modern airway management devices like ET-tube (8 mm I.D. endotracheal tube is recommended), supraglottic airway or Combitube™ etc.

The insertion of the laryngoscope and endotracheal tube is to be made very carefully as if it was a real patient. It is important to grease the tongue part from the oral cavity down to the white stripes simulating the vocal chords using the accompanying water-soluble lubricating gel and lubricate as often as the tongue part feels dry.

Do not pour lubricating gel down the tongue parts but grease the parts with a thin layer of lubricating gel. Furthermore, it is important also to grease the tracheal tube and laryngoscope well with lubricating gel.

Lack of lubrication can destroy the parts, especially the tongue part, but can also make the face piece separate from the tongue part.

Silicone oil or Silicone spray must not be used as the parts can become sticky and almost impossible to clean. The facemask must not be removed by the user. Disassembly should only take place during service or repair.

6.3. Compression ⑦

External chest compression can be practised, and the depth of compression will be displayed on the monitoring instrument in millimetres. It is possible to adjust the stiffness of the chest to simulate a patient with a soft, normal or hard chest.

To ensure that the correct point of compression is applied during chest compression, the instrument also provides indication of incorrect hand positioning, see 7.1.

The carotid pulse can be felt on both sides of the neck during chest compression, or simulated manually by the ECG-control, see 7.2.

6.4. Carotid Pulse

The manikin has an automatic Carotid pulse that will follow the QRS Complex that is set when simulating a pulse using the ECG control in the software application.

Further information about the control can be found in the Directions for use for the Software application.

6.5. Defibrillation and ECG

Caution

Do not interchange the defibrillation electrodes and ECG electrodes on the chest of the manikin.

When using a defibrillator, it is important to follow the manufacturer's instructions for use, particularly the safety rules.

- Before starting the actual training, the instructor should inform the trainees about the proper use of the defibrillator, alerting them to the hazards involved in disregarding the safety rules.
- The manikin is sensitive to moisture. Make sure it does not get wet, as this might cause a short-circuit.

6.5.1. Defibrillation

The manikin allows defibrillation with biphasic and monophasic, manual and -semi-automatic defibrillator or AED's (Automatic External Defibrillators) up to 400 Joule.

6.5.2. Defibrillation electrodes

If standard paddles are used, these should be pressed against the 2 defibrillation electrodes on the manikin's chest. The rhythm can now be read on the defibrillator screen and defibrillation can be performed on the manikin.

Adapter for manual defibrillation is available, which makes the defibrillation electrodes bigger and therefore ensures a better contact between the paddles and the electrodes.

If semi-automatic defibrillators or AED's with adhesive chest electrodes are used, special adapters for the specific defibrillator (simulating the adhesive chest electrodes) should be fitted to the chest defibrillation electrodes on the manikin making it possible to connect the cable from the AED's.

6.5.3. ECG leads

4 ECG connectors are placed on the chest of the manikin making it possible to simulate 3-4 lead ECG and monitor the ECG on a defibrillator and a separate scope.

6.5.4. ECG-Control ⑧

The simulation of the ECG signal can be controlled in the software. Further details can be found in the IFU of the software application.

6.6. I.V. Trainer (optional)

The I.V. Trainer enables practicing placement of an infusion cannula and insertion of fluid (Distilled Water).

Further detailed instructions can be found in the Directions for use for the Ambu I.V. Trainer.

NOTE: During transport the fluid inside the I.V. Trainer must be drained off.

Only use distilled water when filling the I.V. Trainer and when simulating infusion and injection, otherwise the tubes simulating veins and arteries become sticky and the material will decompose.

During storage between training sessions the I.V. Trainer must be covered or placed in the carrying case together with the AmbuMan Advanced. The tubes simulating veins and arteries and the skin are made of natural latex. The arm must therefore not be exposed to sunlight or other U.V. (ultra-violet) rays as the material will become sticky and decompose.

7. Cleaning

7.1. Cleaning of skull, neck and body

The skull, face piece, neck and body can be wiped over with a cloth moistened in a mild detergent and then wiped over once more using a cloth moistened in clean water. When washing, never allow detergent or water to enter the skull, body or instrument unit.

If necessary, cover the connection between body and instrument unit with a cloth.

Marks on the manikin left by lipstick or ball pen can penetrate the material and should therefore be removed as quickly as possible using alcohol.

7.2. Cleaning of clothing

The tracksuit is made of 50% cotton and 50% polyester and is washable at max. 40° C (104° F).

7.3. Cleaning of suitcase

The suitcase can be cleaned using a mild detergent with a cloth or soft brush, rinsed in clean water, and then dried off.

7.4 I. V. Trainer

To clean and maintain the I.V. Trainer, please see the direction for use for the Ambu I.V. Trainer.

8. Quick connection guide

To connect the manikin to a computer device the following steps have to be performed:

1. Turn on the manikin.
2. Connect the WiFi of the Computer device with the network "AmbuW".
3. Open the web browser and type "Ambu.login".

Then follow the instructions on the screen.

If an NFC (near field communication) device is used the connection can be set automatically following the enclosed steps:

1. Turn on NFC on the computer device
2. Place the device on the manikin next to the Power button to connect to the network.
3. Place the device next to the opposite side of the instrument panel to start the software application.

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