

Ambu[®] Man School



Key features

- Ambu Hygienic system to prevent cross contamination
- Compact and portable training manikin
- easy training of chest compressions
- Closed torso
- Easy placement of AED pads
- 10 manikins for easy class training

AmbuMan School

Ambu Man School is a training manikin which offers human appearance dedicated to the training of all procedures for basic life support resuscitation

The manikin can be used for realistic training of external chest compressions and the placement of AED pads in order to get used in the handling of AED's.

The unique hygienic system prevents cross contamination when demonstrating and training of the principals of ventilation through mouth and nose.

The AmbuMan School intended to be used for education of pupils in schools under the supervision of an instructor/ teacher.

Guideline conform

The AmbuMan School allows to perform CPR training following the latest international CPR.

Compact and portable

For an easy transportation & set up possibility the Ambu Man School has a low weight and compact packaging.

Easy AED training

Due to the closed shaped torso, defibrillation pads can be placed and simulated AED training can be performed realistically.

Specifications

Weight	0,8 kg
Length	45 cm
Max. compression depth	6 cm

Material

Straps	Nylon
Base plate	Polycarbonate
Bellow	Thermoplastic Elastomer
Chest	Polypropylene
Head	Silicone
Face piece	Polyvinylchloride
Lung bag	Polyvinylchloride
head bag	Polyethylene

Order information

Item no.	Description
A268 007 000	AmbuMan School (Set of 10)
268 000 571	Face piece AmbuMan School (Set of 10)
268 000 572	Head bag AmbuMan School (Set of 100)
268 000 573	Accessories AmbuMan School (Set of 10)



Closed torso



Placement of AED pads



Ambu hygienic system



Accessories

Ambu BV

Schiphol Boulevard 127
1118 BG Schiphol Airport
The Netherlands
Tel. +31 182 52 60 60
Fax +31 182 52 70 73
www.ambu.nl