

Powerheart® AED G5 Trainer

Our simple-to-use trainer offers a realistic, hands-on rescue experience.

The Powerheart AED G5 Trainer plays a key role in your successful AED programme.

Being Rescue Ready® is more than just having AEDs available, it requires preparation of potential rescuers so they know what to do in a sudden cardiac arrest emergency. The Powerheart G5 Trainer was made to do exactly that. This simple to use trainer has the look and feel of the Powerheart G5. It simulates AED functionality without delivering an electrical shock for a valuable hands-on rescue experience.

- Includes four rescue scenarios
- Instructors can vary AED and rescue conditions to test student responses using the wireless remote
- Scenarios can be paused so instruction can be done throughout the training
- The trainer supports both automatic and semi-automatic AED operation
- Students can practice with reusable adult pads (with or without ICPR) and/or reusable pediatric pads

The G5 AED Trainer provides an easy way to stay in practice.

Owning a G5 Trainer gives you the flexibility to keep skills fresh to ensure your team Rescue Ready®.

Users

- AED/CPR instructors
- Designated workplace responders
- EMS responders
- Anyone wanting to learn how to use an AED

Primary Benefits

Practice. Students gain real world experience to help a sudden cardiac arrest victim. And with CPR feedback, they get a good feel for effective CPR.

Ease of use. Clear, instructive voice prompts guide students through a simulated rescue, and the built-in metronome sets the proper chest compression pace. Configurable for use in a variety of languages for added flexibility.

Preparation. With four pre-programmed rescue scenarios, students will learn what to do in different situations.



Powerheart® AED G5 Trainer



Overview	
Interactive	<ul style="list-style-type: none">· Four pre-programmed rescue scenarios following American Heart Association (AHA) and European Resuscitation Council (ERC) guidelines.· Advance rescue stages and change parameters using an infrared, wireless remote.· Pause scenarios for continuous instruction throughout the training.· Voice and text prompts accurately represent what is heard during an actual rescue situation.· Instructors may issue ICPR feedback corrective prompts for rate, depth and recoil with the wireless remote.
Versatile	<ul style="list-style-type: none">· Change rescue scenarios to introduce different variables· Simulate both shockable and non-shockable rhythms· Monitor student skills during simulations· Includes over 30 language options. Configure the trainer with different languages as needed.
Economical	<ul style="list-style-type: none">· Reusable training electrodes.· Compatible with any type of CPR manikin.· Trainer is powered with three C-cell alkaline batteries; wireless remote uses two AAA alkaline batteries.· The G5 Trainer supports both semi-automatic and fully-automatic AED training.
Safety	<ul style="list-style-type: none">· Cardiac Science AED Trainers will not generate or deliver defibrillation shocks.
Dimensions	<ul style="list-style-type: none">· 8 cm x 32 cm x 24 cm (3.3 in x 12.5 in x 9.5 in)
Part Numbers	<p>190-5020-001: G5 Trainer, Adult Training Pads Package includes: Powerheart G5 Trainer, (1) Set of adult trainer pads, (1) Infra-red wireless remote, (1) User's Guide CD with over 30 language packs, Overlay for automatic version, and Language Configuration Instructions</p> <p>190-5020-002: G5 Trainer, Adult Training Pads with ICPR Package includes: Powerheart G5 Trainer, (1) Set of adult trainer pads with ICPR, (1) Infra-red wireless remote, (1) User's Guide CD with over 30 language packs, Overlay for automatic version, and Language Configuration Instructions</p> <p>XTRPAD004A: Reusable Adult Training Pads</p> <p>XTRPAD006A: Reusable Pediatric Training Pads</p> <p>XTRPAD005A: Reusable Adult Training Pad with ICPR</p>

