

Available in



SmartCase SC1230

Official dealer Benelux: AEDvoordelig.nl

Waterproof
IP54



Active
ventilation



AED
protected
until -20°C

Stainless
Steel



SmartCase SC1230

To support Public Access Defibrillation and provide the best protection for the AED, GT Medicare developed the SC1200 series outdoor AED cabinets. The most effective, straight forward designed AED cabinet, which ensures the **best performance** for the **lowest price**.

The cabinet is equipped with a 24V heating element, active ventilation and a polycarbonate impact-proof window. After operating the door handle the AED can easily be removed and the red LED light will illuminate, accompanied by an intermittent alarm. The SmartCase should be connected to a 230V/ 110V power supply, which the adapter converts to 24V.

Technical Specifications

Type	SC1230
Weight	9 kg
Dimensions	40(h)x40(w)x20(d)cm
Shockproof	✓ IEC61439-5 clause 10.2.101.2/3
Waterproof	✓ IP54
Cooling	✓ Active ventilation
Heating	✓ -20 °C
Corrosion resistant	✓ Stainless steel with coating
Window	✓ Impact-proof polycarbonate
Installation kit	✓ Included (wood & stone)

Power supply

Power 24V Transformer included

Functions

Opening	Door handle/ knob
Light	✓ LED red and blue
Acoustic alarm	✓ +/- 95dB intermittent

Retail Price

€ 425,-

Best performance, most complete,
lowest price...

Amongst others compatible with

